



LIGHT DRIVE 2.1


Designed for all manual wheelchair users, whether occasional, temporary or permanent, LIGHT DRIVE 2.1 turns manual wheelchair into ultralight and high performance powerchair in seconds.


Benoit Systemes has been fitting power add-ons on manual wheelchairs in France since 1994. LIGHT DRIVE 2.1 offers you the possibility to effortlessly add power to your manual wheelchair.



HOW

DOES IT WORK?





Disengage the power add-on in one movement and feel free to self-propel.

Jean-Claude
User since 2016 - diagnosed
with multiple sclerosis in 2013.

My Light Drive helps me all year around to go shopping with my wife, to stroll through shopping centres and, on fine days, to wander around car boot sales and markets. Since I've had my Light Drive, I've become independent. I no longer need somebody with me all the time.





ADVANTAGES OF

LIGHT DRIVE 2.1











LIGHTWEIGHT AND EASILY TRANSPORTABLE

ADAPTS TO ALL WHEELCHAIRS* 'see the compatibility guide

INSTANT INSTALLATION

INDOORS/ OUTDOORS

RANGE UP TO 15 KM



CAN BE TRANSPORTED BY PLANE*



MADE IN FRANCE



REMOVABLE



COLOUR OF YOUR CHOICE

*11.4 Ah battery complies with IATA regulation.



KEY FEATURES



Ergonomic clutch system



Retractable joystick



Anti-tip for safety



Extended colour range



Battery with charge indicator

WIDE RANGE OF COLOUR

METALLICS



CLASSICS

- BLACK
 - K
- MIDNIGHT BLUE
- SKY BLUEPINK

GREEN

TRENDY

BLUERED

KHAKI

WHITE

- GREYSILVER
- CHAMPAGNE GOLD
- ORANGE

LIGHT DRIVE MINI 2.1



A smaller version for children

Light Drive Mini 2.1 can be refitted as the child grows allowing to change wheelchair without changing the power add-on.

Wheelchair equipped with a Light Drive maintains its original manoeuvrability and foldability, it is also easy to transport in any vehicle.

Thanks to Light Drive Mini 2.1, your child will be able to ride easily on all kinds of terrain (steep slopes, grass, gravel).

Reverse gear makes the pedestrian pavement and small step more accessible.



TECHNICAL DETAILS

Our power add-ons can be fitted on most of wheelchairs starting from **32cm** seat width and **40** cm between the wheels.

	Maximum user weight	Power add-on weight
Light Drive 2.1	120 kg	8,8 kg
Light Drive Mini 2.1	120 kg	8,8 kg
Light Drive Plus 2.1	240 kg	16 kg

Max speed ———————————————————————————————————	 10 km/h programmable
Motor —	 2 x 130W motors with brakes
Motor Light Drive Plus 2.1 —————	 4 x 130W motors with brakes
Climbing ability —	— 6° (11%)
Manufacturer's warranty —	24 months for the power add-on
	and the battery

LIGHT DRIVE PLUS 2.1





Light Drive Plus 2.1 has been designed for users of wheelchairs with a weight above 120 kg.

Equipped with four 130 W motors, it will allow you to enjoy the advantages of lightweight power add-ons.

Light Drive 2.1 Plus can be fitted on all foldable wheelchairs. It is provided with a battery of 18Ah and offers a range up to 20 km.



Compatible wheelchairs: a compatibility guide is available to know the models of wheelchairs on which Light Drive 2.1 can be installed, in compliance with general requirement for safety and performance stated in the European Regulation related to medical devices. Find out more



11.4 AH LITHIUM-ION BATTERY COMPLIES WITH APPLICABLE IATA REGULATIONS FOR AIR TRAVEL.

ModelWeightRange- Lithium-ion battery- 3,85 kg- up to 15 km





